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
2- Backup boar fertility prediction software—Ecotext

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3- Intermuscular fat measurement software—Biosoft

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

WHAT parameters of semen analysis? 精液分析有哪些参数?

I want only
“round strawberry flavored
candy”!
我只想要“圆的草莓味的糖果”!

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I want only
“round strawberry
flavored candy”!
我只想要“圆的草莓味的糖果”!

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
I want only
“round strawberry flavored
candy” !
我只想要 “圆的草莓味的糖果” !

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

EK5 Viability & Acrosome integrity 活率和顶体完整性


EK6 Bacterial concentration 细菌浓度

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IMV strenghtening position with easyCyte

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



Bovine牛: $R^2= 0.84$

with high green fluorescence representative of immature cells). A significant relationship ($R^2 = 0.84$, $P < 0.05$) was observed between real and predicted fertility. Once the accuracy of fertility prediction has been confirmed, the model developed in the present study could be used by artificial insemination centers for bull selection or for elimination of poor fertility ejaculates.

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Equine马: $R^2= 0.942$

In conclusion, the results of the present experiment has enabled the determination of a new protocol for the evaluation of stallion semen, combining microscopical observation, computer-assisted motility analysis and flow cytometry, and providing a high level of fertility prediction. Similar investigations are in progress on frozen semen.

(Barrier B. et al. Theriogenology. 2016)


Fig. 3. Relation between the observed field fertility value ("observed fertility") and the predicted fertility value based on the linear regression model ("calculated fertility"). For each stallion ($n = 33$), the observed fertility is presented as a function of the calculated fertility, as well as the minimal and maximal values of the 90% confidence interval ("CI 90% mini" and "CI 90% maxi"). The equation of the regression curve between the observed fertility and the calculated fertility is presented with its r^2 . The adjusted r^2 of the linear regression model that was used to calculate the fertility was 0.942.

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
The newly developed Fertiboar technology is based on the analysis of boar testicular ultrasound pictures which reflect characteristics of the sperm producing tissue (Theriogenology 158 (2020): pp 58-65). A machine learning approach (Convolutional Neural Network) was used to identify regions within the testicular pictures that are key to semen production. Those regions were analysed to calculate parameters describing novel attributes like the homogeneity and echogenicity of the tissue.

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
Source:french alibaba 2017

BioSoft Toolbox®

Biotronics Inc.
Ultrasound Solutions

BioSoft Toolbox® for Swine 4.0.1.1 © 2007-2020 by Biotronics, Inc.

File Tools Options Help

(N) Next Animal

Image Capture Tools
(S) Stream (F) Freeze
(C) Capture Image

Scanner Setup
Depth: 14cm
Probe Size: 13cm

Image Interpretation Tools
☒ (I) IMF
of ROI: 2
ROI Size: 66
☐ w/ Depth Assist
☐ (R) Area
☐ (D) Fat Depth
(A) Assist With Depths

☐ Sort By Scan Time

Animal List
16-4
23-12
24-9
24-10
27-1
28-1
28-6
28-9
28-10
29-1

14-Oct-20 10:57 AM
BioSoft Toolbox®: 4.0.1.0
Biotronics DDI for Exago, a licensed product
Exago Software Version: 01-12

Animal Information
ID 28-10
IMF 0.3
Area 6.83
Fat Depth 0.49
Loin Depth
Weight 232
Sex 2
Breed Berkshire
Farrow Date
Sire
Comments

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BioSoft Toolbox[®] 4.0

Biotronics Inc.
Ultrasound Solutions

B Depth Assist Settings

Measurement Locations

Last Rib - Depth measurement taken anterior of the last rib

☐ Between second to last rib and third to last rib

☐ Between third to last rib and fourth to last rib

☒ From second to last rib to fourth to last rib

10th Rib - Depth measurement taken posterior to the trapezius muscle

☐ Between 10th and 11th Rib

☐ Between 11th and 12th Rib

☐ Between 12th and 13th Rib

☐ Between 11th and 13th Rib


Note: Rib numbers are approximate for market weight pigs and start from the same side as the trapezius muscle.

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BioSoft Toolbox® 4.0

Biotronics Inc.
Ultrasound Solutions

BioSoft Toolbox® for Swine 4.0.1.1 © 2007-2020 by Biotronics, Inc.

File Tools Options Help

(N) Next Animal

Image Capture Tools
(S) Stream (F) Freeze
(C) Capture Image

Scanner Setup
Depth: 14cm
Probe Size: 13cm

Image Interpretation Tools
☐ (I) IMF
of ROI: 2
ROI Size: 66
☐ w/ Depth Assist
☐ (R) Area
☒ (D) Fat Depth
(A) Assist With Depths

☐ Sort By Scan Time

Animal List
16-4
23-12
24-9
24-10
27-1
28-1
28-6
28-9
28-10
29-1

14-Oct-20 10:56 AM
BioSoft Toolbox®: 4.0.1.0
Biotronics DDI for Exago, a licensed product
Exago Software Version: 01-12

Animal Information
ID 28-10
IMF 0.3
Area 6.83
Fat Depth 0.50
Loin Depth 2.54
Weight 232
Sex 2
Breed Berkshire
Farrow Date
Sire
Comments

Images Image Data In English Units
1 2 3 4 5

Animals: 35 Animals w/ Images: 35 Pixel Value: 0 Hardlock Count: 6144 File Version: 4.0.1.0 Session Folder: Z:\BioQProc\Previous Years\2020\AJC Berkshire...

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BioSoft Toolbox® 4.0

Biotronics Inc.
Ultrasound Solutions

BioSoft Toolbox® for Swine 4.0.1.1 © 2007-2020 by Biotronics, Inc.

File Tools Options Help

(N) Next Animal

Image Capture Tools

(S) Stream (F) Freeze

(C) Capture Image

Scanner Setup

Depth: 14cm

Probe Size: 13cm

Image Interpretation Tools

☐ (I) IMF

of ROI: 2

ROI Size: 66

☐ w/ Depth Assist

☒ (R) Area

☐ (D) Fat Depth

(A) Assist With Depths

☐ Sort By Scan Time

Animal List

- 16-4
- 23-12
- 24-9
- 24-10
- 27-1
- 28-1
- 28-6
- 28-9
- 28-10**
- 29-1

14-Oct-20 10:56 AM

BioSoft Toolbox®: 4.0.1.0
Biotronics DDI for Exago, a licensed product
Exago Software Version: 01-12

Exago

Animal Information

ID	28-10
IMF	0.3
Area	7.16
Fat Depth	0.50
Loin Depth	
Weight	232
Sex	2
Breed	Berkshire
Farrow Date	
Sire	
Comments	

Images Image Data In English Units

1 2 3 4 5

Animals: 35 Animals w/ Images: 35 Pixel Value: 0 Hardlock Count: 6145 File Version: 4.0.1.0 Session Folder: Z:\BioQProc\Previous Years\2020\AJC Berkshire...

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IMV TECHNOLOGIES. FOR LIFE.

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Thanks for your attention!!

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Copyright © 2007-2020 by Biotronics, Inc.